Medtronic

Sophono™

Magnetic Implant

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Implanto magnetico

Magnetimplantat

Implante magnético

Magnetisch implantaat

Magnetimplantat

Magneetti-implantti

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Implante magnético

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EN

Indications for use

The Sophono Magnetic Implant is intended for use with the Processor (patients 5 years of age and up) for the following patients and indications:

- Patients with conductive or mixed hearing losses, who can still benefit from amplification of sound. The pure tone average (PTA) bone conduction (BC) threshold for the indicated ear should be better than 45 dB HL (measured at 0.5, 1, 2, and 3 kHz).
- Bilateral fitting is applicable for most patients having a symmetrically conductive or mixed hearing loss. The difference between the left and right side BC thresholds should be less than 10 dB on average measured at 0.5, 1, 2, and 4 kHz, or less than 15 dB at individual frequencies.
- Patients who have a profound sensorineural hearing loss in one ear and normal hearing in the opposite ear, who for some reason will not or cannot use an AC CROS. The pure tone average (PTA) air conduction (AC) threshold of the hearing ear should be better than 20 dB HL (measured at 0.5, 1, 2 and 3 kHz).

Device description

The Sophono Magnetic Bone Conduction Hearing System is a family of Processors and accessories that operate on the principle of bone conduction of sound vibrations.

The Sophono Magnetic Implant is configured to be magnetically attracted to the Processor. The Magnetic Implant holds the Processor against the head, and vibration is transduced through direct contact with the patient's skin and the bone below.

The Sophono System is designed for use for those patients with conductive hearing loss, those patients who have sensorineural hearing loss up to 45 dB in combination with their conductive loss, and single sided deafness as defined in the indications for use. The prescriptive formula and adjustments available to the audiologist in the software allow for programming the Sophono System for individual patient hearing loss.

Contraindications

Contraindications for the Sophono Magnetic Implant connected to the Processor include factors that would cause a clinician to refer the patient for medical assessment which will temporarily, or in some cases permanently, halt the process of Processor fitting. These factors are:

- A hearing loss of sudden onset;
- A rapidly progressing hearing loss;
- Pain in either ear;
- Tinnitus of sudden recent onset, or unilateral tinnitus;
- Unilateral or markedly asymmetrical hearing loss of unknown origin;
- Vertigo (e.g. dizziness)

Warnings

The following section contains a description of all physician related Warnings associated with the Sophono Magnetic Bone Conduction Hearing Systems. Additional patient and use related Warnings can be found in the **Sophono Bone Conduction Hearing System Instructions For Use**. Knowledge of the following information is essential for safe operation of the device. Specific Warnings pertaining to physician instructions are included.

- Thoroughly read the **Sophono Magnetic Implant Instructions For Use** and the **Sophono Bone Conduction Hearing System Instructions For Use** before prescribing and using the Sophono Bone Conduction Hearing System.
- To avoid potential performance and / or safety issues, adhere carefully to patient selection criteria.
- Blunt force trauma in the area of the Processor may result in personal injury and damage to the magnetic implant or Processor.
- Do not perform Magnetic Resonance Imaging (MRI) on a patient wearing a Processor, Attract™ Magnetic Spacer, Headband, or Softband. Severe injury may occur in an MRI scan room or when undergoing an MRI scan with the external components in place. External components of the device must be removed prior to scanning.
- Implanted subjects may safely receive an MRI scan with the Sophono (M) Magnetic Implant only under very specific conditions. Severe injury can occur if these conditions are not followed. Consult with a physician prior to undergoing an MRI and refer to information in this IFU.
- Non-clinical testing demonstrated the Sophono (M) Magnetic Implant is MR Conditional and can be scanned safely using MRI only under the following conditions:
 - remove all external components including the Processor, Attract™ Magnetic Spacer, Headband or Softband before entering the MR environment
 - static magnetic field of 3 Tesla or less
 - spatial gradient field of 720 Gauss/cm or less
 - maximum whole-body averaged specific absorption rate (SAR) of 4 W/kg in the First Level Controlled Mode for a maximum scan time of 15 minutes of continuous scanning.

- After MRI exposure the Attract™ Magnetic Spacer and Strength should be re-evaluated and adjusted per the IFU if needed.
- The Sophono Magnetic Implant and packing are provided sterile and are intended for single use only. To reduce the potential for infection, do not re-sterilize, do not re-use. Inspect the packaging to be sure that it is intact and undamaged prior to use. To reduce the potential for infection, do not use the product if the pouch is opened or damaged.
- Purchased screws may be provided non-sterile and if so must be sterilized prior to implantation according to the manufacturer's instructions. Failure to sterilize appropriately can result in serious infection.
- Both the Alpha (S) and Alpha (M) configurations of the device contain magnets which could impact the performance of magnetically sensitive medical devices such as a pacemaker, an implantable cardioverter defibrillator (ICD), or ventricular shunt. The magnets of the softband/headband of the Alpha (S) and the Attract™ Magnetic Spacer of the Alpha (M) should therefore be kept away from such devices.

Precautions

The following section contains a description of all physician related Precautions associated with the Sophono Magnetic Bone Conduction Hearing Systems. Additional patient and use related Precautions can be found in the Processor Instructions For Use. Knowledge of the following information is essential for safe operation of the device. Specific Precautions pertaining to physician instructions are included.

- The Sophono Magnetic Implant is intended for single-patient use only. Do not use any of these products if they have previously been placed in or with another patient. Reuse of the magnetic implant carries a serious risk of severe infection and other major medical complications.
- Patients with a pacemaker or ICD should consult their care provider for safe use of the Sophono system with their cardiac medical device.
- Keep Sophono System components away from magnetic data storage media and electronic devices. The Processor, Magnetic Implant, Headband, and Softband all contain magnets that can potentially damage magnetic data storage media and electronic devices.
- To avoid potential of damage to the device, do not store the Sophono System near magnets or magnetic devices other than the Attract™ Magnetic Spacer and implant.
- Refer to the original manufacturer's instructions for use for specific directions on how to use the implantable screws.

Pre-surgical training requirements

Before attempting to implant the Sophono Magnetic Implant, surgeons are required to complete training provided by Medtronic. During the training they will receive guidance regarding the procedure, a thorough understanding of the device and surgical procedure and other important information. Contact Medtronic Customer Service to arrange training.

Surgical risks

The implantation of the Sophono Implant requires a surgical procedure. As with any surgical procedure involving a general anesthetic, certain risks exist. Those risks include the risk of infection from poor hygiene or non-sterile technique, irritation of tissues during the implantation procedure, unintended cuts and abrasions, and irritation of tissues due to unintended exposure to implant materials. The implant is made of materials considered to be biocompatible but in the case of unintended damage to the device, a patient could come in contact with substances that could generate an adverse reaction.

Patient selection criteria

- Patients to be fitted with the Processor (S) and its Headband or Softband may be any age.
- Patients to be fitted with the Processor (M) and its Magnetic Implant should be 5 years old or older.
- Patients and/or patient guardians must be able to understand that the Sophono System will not restore hearing to normal and will not prevent additional hearing impairment resulting from conditions such as disease, old age hearing loss, genetic hearing loss, ototoxic drugs or noise induced hearing loss.

Possible adverse events

It is possible that a reaction to the materials of the Magnetic Implant could occur. Such a circumstance may necessitate that use of the device be discontinued.

Risks that may occur during surgery

The risks involved in surgery and complications related to surgery include the normal risks of the operative procedure and general anesthesia. These include pain, blood clots and/or circulatory problems, infection, cardiac arrest, and even death.

Specific risks related to operating on the mastoid include infection, inflammation, numbness or stiffness about the ear, dural exposure, Cerebral Spinal Fluid (CSF) leak, perforation of sigmoid sinus, subdural hematoma, disturbance of taste or balance, and noticeable changes in head noise.

Risks that may occur after surgery

Potential risks after surgery include inflammation and complications associated with general anesthesia. Potential risks related to the implant following surgery include infection, reparative granuloma and inflammatory reaction.

Specific risks related to operating on the mastoid include infection, inflammation, numbness or stiffness about the ear, dural exposure, Cerebral Spinal Fluid (CSF) leak, perforation of sigmoid sinus, subdural hematoma, disturbance of taste or balance, and noticeable changes in head noise.

It is possible that a reaction to the implant materials and or failure of the implanted device can occur and may necessitate an operation to remove and/or replace the implant.

Preparing for surgery

Pre-Surgical antibiotic prophylaxis

Antibiotic prophylaxis should be considered if the surgeon determines that this surgical procedure may be associated with a considerable risk of infection or if a post-operative infection would cause a serious hazard to the patient's recovery or well being.

Marking the incision

Prior to entering the sterile field, use the non-sterile template and a surgical marking pen to mark the preferred position of the incision on the patient's scalp to avoid placing incision over the top of the implant location. Refer to Figure 1, for orientation and position.

Preparing the surgical site

The hair should be clipped and the scalp shaved at the incision site immediately before surgery. With the incision site skin free of gross contamination, apply antiseptic (such as iodophors, alcohol-containing products, or chlorhexidine gluconate) in the area of the proposed incision to an area large enough to allow for the extension of the incision or creation of a new incision or drain sites (if necessary).

Sophono (M) system compatible screws

The Sophono (M) implant requires 5, self-drilling/self-tapping, craniomaxillofacial, titanium screws to be used to affix the implant to the skull. The implant bracket holes are sized to accept commercially available craniomaxillofacial titanium screws. Maximum screw diameter is less than 2.0 mm and the minimum screw head diameter is greater than 2.45 mm. Screws must be 4 mm in length (3 mm for emergency screws). These screws are NOT included as part of the Sophono (M) implant system and must be obtained separately.

Preferred Craniomaxillofacial Titanium Screw sizes	1.6 to 1.7 mm diameter x 4.0 mm Self tapping	1.6 to 1.7 mm diameter x 4.0 mm Self drilling	1.8 to 1.9 mm diameter x 3.0 mm Emergency
---	---	--	--

If the initial attempt to place a screw does not work and the screw hole becomes "stripped" an emergency "rescue" screw can be threaded into the stripped hole. The largest size screw that can be used is 1.9 mm in diameter.

The screws must be sterile prior to implantation. If the screws are not purchased sterile, refer to the screw manufacturer's instructions for sterilization.

Standard bone bed implant procedure

Making the incision

The proper position for the implant is 6 cm from the ear canal at a 45-degree angle posterior and superior to the ear canal. The incision should be made 7.5 to 8 cm from the external ear canal (1.5 to 2 cm from the location of the implant so the incision is not directly over the implant). The Surgical Template will aid the surgeon in ensuring proper placement by marking the skin with a surgical marking pen. After making the incision, expose the skull where the implant will be placed. The incision must be large enough to allow placement of the implant. See Figure 1.

Creating a bone bed (optional)

The implant is approximately 2.6 mm thick and optimal tissue thickness over the implant is 4 mm to 5 mm. The Surgical Template is provided with 6 mm and 4 mm gauges for determining the appropriate skin flap thickness. Bone beds are recommended for the implant if the tissue thickness over the implant is less than 3 mm. See Figure 2.

To create a bone bed, mark the position of the implant magnets on the bone using the Surgical Template and drill the bone bed to a depth of approximately 3 mm. The implant magnets have a diameter of 10 mm. The Implant Template is useful in ensuring that the bone beds are of sufficient depth and diameter.

Placing the implant

After creation of the bone bed, install the magnetic implant using the specified screws such that the magnets of the implant are fully seated in the bone beds. See Figure 2.

Place the center screw first. Place one screw in each of the bracket holes to securely attach the implant to the skull. The surface of the implant should be nearly flush with the surface of the mastoid. See Figure 3.

Checking the skin thickness and closing the incision

If the tissue flap is thicker than 6 mm, thin down a 10 mm wide strip of tissue to 5 mm directly over the implant magnets. Suture the incision closed.

Minimally invasive implant procedure

Making the incision

The proper position for the implant is 6 cm from the ear canal at a 45-degree angle posterior and superior to the ear canal. The incision should be made 7.5 to 8 cm from the external ear canal (1.5 to 2 cm from the location of the implant so the incision is not directly over the implant). The Surgical Template will aid the surgeon in ensuring proper placement by marking the skin with a surgical marking pen. After making the incision, expose the skull where the implant will be placed. The incision must be large enough to allow placement of the implant. See Figure 1.

Creating the skin flap

Clean bone under skin flap (use periosteal elevator). Preserve periosteum where possible (for healing and protection). Assure implant template fits well on exposed bone; not under the incision line.

Placing the implant

Place the flat side of the implant against the surface of the skull. See Figure 4. Place center screw first. Gently depress implant arms so implant is flush with bone surface. Place and secure other screws. Make sure implant is flush and all screws secure.

Checking the skin thickness and close the incision

Using template gauge, check for skin thickness of 4-6 mm. If the tissue flap is thicker than 6 mm, thin down a 10 mm wide strip of tissue to 5 mm directly over the implant magnets. Move template gauge sideways to make sure skin is thinner than gauge opening. If necessary thin soft tissue evenly, (do not create a divot over the implant). Suture periosteum and skin.

Placing the Processor

Approximately 4 weeks after surgery, after the incision site has healed and all signs of inflammation are gone, the patient may be fitted with the Processor. Refer to the **Sophono Bone Conduction Hearing System Instructions For Use** for information.

Removing the magnetic implant

Making the incision

Locate the position of the implant. The implant is approximately 6 cm from the ear canal at a 45-degree angle posterior and superior to the ear canal. The incision should be made 7.5 to 8 cm from the external ear canal (1.5 to 2 cm from the location of the implant so the incision is not directly over the implant). After making the incision, expose the skull where the implant has been placed. The incision must be large enough to allow removal of the implant. See Figure 1.

Removing the implant

Remove each of the 5 screws holding the magnetic implant in place. It may be necessary to excise bone or tissue growth occurring over the face of the implant. Once the screws are removed the implant can be lifted carefully from the bone or bone beds. Suture incision closed.

Magnetic fields

The Sophono (M) Magnetic Implant creates a permanent magnetic field. A permanent magnetic field is not the same as, and should not be compared to, an electromagnetic field such as that of a mobile telephone or a high-voltage transmission line.

MRI safety information

Non-clinical testing has demonstrated the Sophono (M) Magnet Implant is MR Conditional. Patients who have the Sophono Magnetic Implant can undergo an MRI scan under very specific conditions in both 1.5 and 3-tesla magnetic fields with the magnetic implant intact. Consult a physician prior to undergoing an MRI. The patient can be scanned safely using MRI under the following conditions:

- remove all external components including the Processor, Attract™ Magnetic Spacer, Headband or Softband before entering the MR environment
- static magnetic field of 3 Tesla or less
- spatial gradient field of 720 Gauss/cm or less
- maximum whole-body averaged specific absorption rate (SAR) of 4 W/kg in the First Level Controlled Mode for a maximum scan time of 15 minutes of continuous scanning.

In non-clinical testing the Alpha (M) Magnetic Implant produced a maximum temperature rise less than 3.2°C during 15 minutes of continuous MR scanning in the First Level Controlled Mode at a maximum whole-body averaged SAR of 4 W/kg.

The computed implant temperature increase in response to the worst-case, time-averaged gradient field (94.7 Tesla/sec) for the 15 minute possible exposure for a series of clinical MRI scans is less than 2.6°C.

MR image artifact extends approximately 5 cm relative to the size and shape of the magnetic implant in a gradient echo scan at 3 Tesla.

After MRI exposure the Attract™ Magnetic Spacer Strength should be re-evaluated and adjusted per the IFU if needed.

MRI patient positioning

To minimize demagnetization of the internal magnets, it is recommended that the patient's head be positioned along the long axis of the MRI system and that the patient is instructed to avoid head movement while on the patient table. If the head is tilted along the long axis of the patient table during a scan, the internal magnets may be demagnetized and a stronger external Attract™ Magnetic Spacer may be needed to keep the Processor in place. See Figures 5, 6, and 7.

To minimize demagnetization avoid facing the side of the patient's head (ear) toward or away from the bore of the MRI system at any time. As much as possible, maintain 90° between the face of the magnetic implant and the static magnetic field of the MRI system. See Figures 5, 6, and 7.

Do not undergo scanning with Magnetic Resonance Imaging (MRI) while wearing a Processor, Attract™ Magnetic Spacer, Headband or Softband. Severe injury may occur in an MRI scan room or when undergoing an MRI scan with the external components in place. External components of the device must be removed prior to scanning.

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Medtronic B.V.
Medtronic Australasia Pty Ltd
Medtronic of Canada Ltd